

# Tecotherm Neo

Thermoregulation and Monitoring



Total Body Cooling for infants  
suffering from Hypoxic Ischaemic  
Encephalopathy (HIE)

# Tecotherm Neo

## Thermoregulation and Monitoring


Total Body Cooling is a therapy for infants suffering from Hypoxic Ischaemic Encephalopathy (HIE) after Birth Asphyxia. The original Tecotherm was used in the biggest Randomised Control Trial for cooling infants with HIE, which showed improved neurologic outcomes for the cooled group<sup>1</sup>. Now building on from that success, the Tecotherm Neo takes Total Body Cooling to the next level.

The innovative servo controlled design has instant feedback monitoring the infant's temperature every 2 seconds and making minute changes to the fluid to ensure that the infant's temperature remains stable. Simple alarms are on hand should the temperature deviate more than 0.5°C from the set temperature. The large colour display shows the temperature clearly and colour changes alert nursing staff to any problem.

The Tecotherm Neo can be programmed for a complete treatment cycle, the user can set the target temperatures, duration and even the rate of re-warming / cooling making it easy and simple to set up and use. Changes can be made at any time and all changes, set and measured parameters are recorded on to a memory card every minute for later analysis.



Stand by screen is easy to see across the room.



Clinician can set the duration of treatment, the rate of cooling / warming as well as the target temperature.



A lightweight and portable device, the Tecotherm Neo is the gold standard in Total Body Cooling.

## Modes of Operation

The Tecotherm Neo has 3 modes of operation:

### **Servo-Control Complete Treatment Mode**

Clinician can set the duration of treatment, the rate of cooling / warming as well as the target temperature. The Tecotherm Neo will complete the whole cycle without interruption and then maintain the final temperature.

### **Servo-Control Mode (Constant Rectal Temperature)**

Clinician can set the target temperature, the time to get to target temperature and the length of time to maintain the target temperature. The Tecotherm Neo will complete the whole cycle without interruption.

### **Constant Mattress Temperature Mode**

No Servo Control, useful for research work and where control of the mattress needs to be stable, or where experienced staff are available to monitor the patient.

## The Tecotherm Neo takes Total Body Cooling to the next level

### Flexible

- Clinician can set and adjust parameters at the start and throughout the treatment
- Customise and store treatment profiles (up to nine can be stored)
- Compatible with reusable and disposable accessories
- Optional use of skin probes

### Easy to Use


- Large colour screen is easy to see across the room
- Small lightweight portable device
- The device restarts from the same settings in the event of power failure
- Maintains normothermia after rewarming

### Compatible


- Data storage for allowing simple analysis of temperature profiles
- Export individual patient records to USB



## Example of patient data collected from Tecotherm Neo



# Technical Description and Specifications

<b>Options</b>	Cooling & Warming	<b>Ambient conditions</b>	
<b>Patient weight max.</b>	<50kg	Operation / Treatment Ambient Temperatures	+ 5°C to + 27°C
<b>Dimensions</b>	375mm (width) x 190mm / 215mm (height) x 310mm (depth)	Operation / Treatment Relative Humidity	10% to 75%, not condensating
<b>Weight without accessories</b>	approx. 7.2kg	<b>System safety</b>	
<b>Central cooling module</b>	Thermoelectrically based module	Protection class	Class 1, Risk Class II b, Type BF 
<b>Treatment temperature control ranges</b>	Mattress + 12°C to + 39°C Rectal BCT + 30°C to + 38°C	Standards	DIN EN 60601-1 DIN EN 60601-1-2 DI 60601-1-6 DIN EN 60601-1-8 DIN EN 60601-1-10 DIN EN 60601-2-35 E/F
<b>Hydraulic circulation system</b>			
System pressure max.	0,5 bar		
Flow rate without / with mattress	500 ml / min (shorted) / up to 300 in use		
Internal fluid reservoir capacity	approx. 250 ml		
Circulating fluid	Sterile Water		
Connectors / Couplings	Quick Disconnect Couplings		
Fill up / Refill	Fill up set		
<b>Electrical parameters</b>			
Supply voltage / Mains	100-130V and 200-240V, 50-60 Hz		
Power consumption	max. 350 W		
<b>Patient safety / alarms</b>			
Lower temperature alarm limit	+10°C		
Upper temperature alarm limit	+41°C		
Set temperatures, lower limit	+12°C		
Set temperatures, upper limit	+39°C		

## References

1. Azzopardi et al (2009) Moderate Hypothermia to Treat Perinatal Asphyxial Encephalopathy, N Engl J Med 2009; 361:1349-1358

The details given in this brochure are correct at time of going to press. The company reserves the right to improve the products shown.

500-126(1) 2018-11



Distributed by

### Inspiration Healthcare Limited

2 Satellite Business Village, Fleming Way,  
Crawley, West Sussex, RH10 9NE  
United Kingdom

**T** +44 (0)1455 840 555

**F** +44 (0)1455 841 464

**E** info@inspiration-healthcare.com

**W** inspiration-healthcare.com