

# Resuscitation Unit 104 E type

*Prompt and safe resuscitation is possible.*

The Resuscitation Unit ensures safe resuscitation.



## Specifications

Atom's code. 60982 Resuscitation Unit 104 E type (Stand-alone type)	
General	Pressure indication: -2~8kPa (-20~80cmH <sub>2</sub> O)
Suction unit	Suction pressure control range: 0~-26.7kPa (0~-200mmHg)
Oxygen blender / flow meter unit	Oxygen concentration range: 21~100% Flow rate range: 0~15L/min
Resusci Flow unit	Maximum pressure limit (P <sub>MAX</sub> ) setting range: 2.0~5.9kPa (20~60cmH <sub>2</sub> O) Default setting: 3.9kPa (40cmH <sub>2</sub> O)
Accessories	Patient Circuit .....1    Test Bag, 50mL .....1    Piping Connecting Hose for Air ...1    Piping Connecting Hose for Oxygen ...1 Outlet Connecting Hose for Air/Oxygen ...1    Oxygen Supply Hose .....1    Face Mask for Neonatal .....1 Face Mask for Infant .....1    Oxygen Humidifier Jar Holder .....1    Humidifier Bottle OX-230 .....1 Outlet Connecting Hose for Suction .....1    Suction Jar Holder .....1    Suction Canister, 500mL .....1 Dispo Liner, 500mL ....1    Suction Hose ....1    Silencer (for E-type) .....1    Operation Manual .....1

〈Note〉 T-Piece Resuscitation Unit with Blender and Suction Unit of ejector type to make suction operation possible where no suction piping is available.

## Option & Supplies

Atom's code. 61616 **MF-Rail** For F-Rail (Left & Right)

Atom's code. 61619 **Suction Canister, 500mL** For S401D


Atom's code. 60355 **Dispo Liner, 500mL**  
For S401D, Suction Canister (30pcs/box)

① Atom's code. 60386 **Patient Circuit** For Resusci Flow (5pcs/box)

② Atom's code. 60298 **Test Bag, 50mL** (5pcs/pkg)

③ Atom's code. 60145 **Face Mask** For Infant (20pcs/box)

④ Atom's code. 60146 **Face Mask** For Neonatal (20pcs/box)



Connect Medizintechnik GmbH  
Gspanngasse 4 | A-2130 Mistelbach  
Tel: +43 2572 32400 | Fax: +43 2572 20404  
office@connect-medizintechnik.at  
www.connect-medizintechnik.at

**Manufacturer: ATOM MEDICAL CORPORATION**  
3-18-15, Hongo, Bunkyo-ku, Tokyo, JAPAN  
<http://www.atomed.co.jp>

**Exporter: ATOM MEDICAL INTERNATIONAL, INC.**  
Iwakata Bldg., 3-18-16, Hongo, Bunkyo-ku, Tokyo, JAPAN  
Phone: + 81 3 3815 2941 Fax: + 81 3 3812 9670/4080